




**TAMPA, FLORIDA**  
DECEMBER 5, 2020

# 3-Year Results from the MIMICS-2 Study

# DISCLOSURES

# BioMimics 3D<sup>®</sup>: The Swirling Flow<sup>®</sup> Stent



- Helical centerline
- Simple, accurate placement using standard delivery system




- Imparts non-planar curvature to stented femoropopliteal segment<sup>1</sup>
- Improved biomechanical performance compared to straight stents<sup>1</sup>
- 0% stent fracture MIMICS-2 IDE Study

*BioMimics 3D and Swirling Flow are Registered Trademarks of Veryan Medical Limited*

The BioMimics 3D Vascular Stent System has FDA, PMDA and CE Mark approval. Not available for sale in Japan  
CAUTION: Federal law restricts this device to sale by or on the order of a physician.

# BioMimics 3D: Helical Centerline Promotes Swirling Flow



- Swirling flow increases wall shear stress (WSS) on endothelial cells
- WSS naturally protects against atherosclerosis and intimal thickening <sup>(2)</sup>
- Increased WSS has been shown <sup>(3)</sup> to provide an antiproliferative effect after stenting, without the need for a drug



1. Malek AM et al, JAMA 1999;252:2035–2042

2. Caro CG, Arterioscler Thromb Vasc Biol 2009, 29:158-161

3. Caro CG et al, J R Soc Interface 10: 20130578

# In Vivo CFD Modelling of Swirling Flow in the Stented Segment


LifeStent (left) and BioMimics 3D Stent (right)




# Improving the Biomechanics

## Bent knee cadaver study


### BioMimics 3D



### Straight Stent




# MIMICS Clinical Programme



# MIMICS Randomized Controlled Trial

<b>8 Investigational Sites</b> <i>Corelab (DUS; angiography; Xray)</i>		
FIM Lead-in	N=10 BioMimics 3D	
Prospective Randomization	BioMimics 3D N=50	Straight nitinol stent N=26
24-Month Primary Patency (p=0.05)	<b>72%</b>	<b>55%</b>
Freedom from CDTLR 12-24 months (p=0.03)	<b>91%</b>	<b>76%</b>

## Freedom from CDTLR Landmark Analysis<sup>2</sup>




**Provided the first clinical proof supporting the durable outcome benefit arising from the BioMimics 3D stent compared to a straight nitinol stent<sup>1, 2</sup>**

1. Zeller T et al; *Circ Cardiovasc Interv.* 2016

2. Sullivan TM et al; *Int J Vasc Med.* 2018



# MIMICS Clinical Programme



# Evaluation of Safety and Effectiveness of the BioMimics 3D Stent System in the Femoropopliteal Arteries of Patients with Symptomatic Peripheral Arterial Disease

MIMICS 2

## • Primary Endpoints

- Safety: composite of death, major amputation or CDTLR through 30 days
- Effectiveness: primary patency at 12-months

## • Follow-up: 3 years

## • 43 investigational sites enrolled 271 subjects

- US: 31 sites N = 162
- Germany: 6 sites N = 78
- Japan: 6 sites N = 31

## • Study Principal Investigators

- Timothy M. Sullivan, MD Minneapolis, MN, USA
- Thomas Zeller, MD Bad Krozingen, Germany
- Masato Nakamura, MD Tokyo, Japan

## Independent Core labs:

- ultrasound; angiography; X-ray

## Independent Clinical Event Committee

- adjudication

# Baseline Patient Demographics

		N= 271 Subjects
<b>Age</b>	Mean years $\pm$ SD (N)	68.4 $\pm$ 9.5 (271/271)
<b>Gender</b>	Male / Female	180 (66.4%) / 91 (33.6%)
<b>Risk Factors</b>	Diabetes Mellitus	45.4% (123/271)
	Hypertension	90.0% (244/271)
	Hypercholesterolemia	81.9% (222/271)
	Smoker Current / Former	80.8% (219/271)
<b>Coronary Revascularization</b>	Previous Percutaneous or Surgical	43.2% (117/271)
<b>Previous Peripheral Intervention</b>	None in target vessel	98.2% (266/271)
<b>Rutherford Category</b>	1	0% (0/271)
	2	26.9% (73/271)
	3	67.5% (183/271)
	4	5.2% (14/271)
	5	0.4% (1/271)
<b>Ankle Brachial Index</b>	Mean $\pm$ SD (N)	0.70 $\pm$ 0.20 (257/271)

# Baseline Angiography and QVA

Core Laboratory Data		N= 271 Subjects
Reference Vessel Diameter (mm)	Mean ± SD	5.2 ± 0.9 (269/271)
Lesion Type <sup>1</sup>	De novo	100% (271/271)
Lesion Location in Femoropopliteal Artery	Prox	11.5% (31/270)
	Mid	48.1% (130/270)
	Distal	40.4% (109/270)
Diameter Stenosis (%)	Mean ± SD	77.8 ± 18.3 (269/271)
Lesion Length (mm)	Mean ± SD	81.2 ± 38.4 (269/271)
Total Occlusion (%)		30.0 (81/270)
Calcification (%)	None - Mild	54.1 (146/270)
	Moderate - Severe	45.9 (124/270)
Run-off (%) - 1 or more patent tibial artery (<50% stenosis)		98.8 (237/240)

<sup>1</sup> Investigator-reported

		<b>N= 271 Subjects</b>
<b>BioMimics 3D stents placed<sup>1</sup></b>	<b># Stents / N</b>	<b>305 / 271</b>
	# Subjects with 1 stent	87.5% (237/271)
	# Subjects with 2 stents	12.5% (34/271)
<b>Stented Segment Length<sup>2</sup></b>	<b>Mean ± SD (mm)</b>	<b>112.3 ± 36.3 (269/271)</b>
<b>Diameter Stenosis<sup>2</sup></b>	Pre-stent % ± SD	77.8 ± 18.3 (269/271)
	Post-stent % ± SD	12.6 ± 7.5 (269/271)
<b>Dissections<sup>2</sup></b>	No Dissection	97.8% (263/269)
	Type A-C	2.2% (6/269)
	Type D-F	0% (0/269)
<b>Device Success</b>		<b>100% (271/271)</b>
<b>Technical Success</b>		<b>100% (269/269)</b>

Technical Success: Core Lab determined ≤50% residual diameter stenosis (in-stent) at end of index procedure

Device Success: Successful delivery of System; placement of stent and retrieval of System

1 Investigator-reported

2 CoreLab-reported

# Primary Endpoint: Safety

Composite of CEC-adjudicated Major Adverse Events through 30 days, including death, any major amputation performed on the target limb, or Clinically-Driven Target Lesion Revascularisation

	Performance Goal	Rate (n/N) [95% CI]
Freedom from MAE through 30 days	>88%	<b>99.6%</b> (268/269) [97.7%, 100%]
<b>Primary safety endpoint</b>		<b>Achieved</b>

# Primary Endpoint: Effectiveness


Primary stent patency rate at 12 months.

	Performance Goal	Rate (n/N) [95% CI]
Primary stent patency	>66%	<b>73.1%</b> (182/249) [67.3%, 78.2%]
<b>Primary effectiveness endpoint</b>		<b>Achieved</b>

Patency was defined as no significant reduction in luminal diameter (< 50% diameter stenosis) since the index procedure.

Loss of patency was determined by an independent core laboratory when the peak systolic velocity ratio (PSVR) exceeds 2.0, or where angiography revealed > 50% diameter stenosis, or where the subject had a CDTLR.

# Freedom from Clinically-Driven TLR at 3 Years = 81%



**1-year = 89%**

**2-year = 84%**


**3-year = 81%**

*\*Core Lab adjudicated, clinically-driven TLR with objective evidence  
Subjects are censored at their last known follow-up, or at time of study exit (withdrawal or lost to follow-up) or death*



# Freedom from Clinically-Driven TLR at 3 Years Lesion Calcification

Calcification % at baseline	None - Mild	54.1 (146/270)
	Moderate - Severe	45.9 (124/270)



None/Mild = 78%  
Mod/Severe = 82%

**No statistically significant difference between the groups**

*Subjects are censored at their last DUS follow-up, or at time of study exit (withdrawal or lost to follow-up) or death.*


*Subjects who are PATENT at DUS follow-up are censored at the end of follow-up window*

*Data on file at Veryan Medical*

# Freedom from Clinically-Driven TLR at 3 Years CTO vs. Stenosis

Total Occlusion % at baseline

30.0 (81/270)




CTO = 76%  
Stenosis = 82%

**No statistically significant difference between the groups**

Subjects are censored at their last DUS follow-up, or at time of study exit (withdrawal or lost to follow-up) or death.  
Subjects who are PATENT at DUS follow-up are censored at the end of follow-up window

# Freedom from Clinically-Driven TLR at 3 Years Tercile of Lesion Length

Lesion Length (mm)

Mean  $\pm$  SD81.2  $\pm$  38.4 (269/271)

Tercile 1 = 82%

Tercile 2 = 81%


Tercile 3 = 77%

**No statistically significant difference between the groups**

Subjects are censored at their last DUS follow-up, or at time of study exit (withdrawal or lost to follow-up) or death.

Subjects who are PATENT at DUS follow-up are censored at the end of follow-up window

# Mimetic/DES Study Comparisons-Freedom from CDTLR



Comparable outcomes to DES and Supera despite more challenging lesions and without the need for lesion preparation

# MIMICS-2 Results

- Reproducible, rigorous, high quality data from US, Japan and Europe
- 81% freedom from CDTLR at 3 years
- Comparable outcomes to DCB, DES and Supera despite more challenging lesions and without the need for lesion preparation – providing ease-of-use simplicity and long-term benefits
- 0% stent fracture

## BioMimics 3D



MIMICS-2 study shows continuing benefit of BioMimics 3D at **3 Years**, even in challenging cases

# Conclusion

BioMimics 3D  
reduces the burden of  
re-intervention for the  
patient and the health  
care system<sup>1, 2</sup>



# Post Approval Study

A prospective, multicenter observational study evaluating the safety, effectiveness and device performance of the BioMimics 3D Vascular Stent System within a real-world clinical population of patients undergoing femoropopliteal intervention.

- 40 sites across US
- 500 total patients
- CEC and Independent Core Lab
- First enrollments targeted for Q4 2020
- Coordinating Investigators:
  - Dr. Sahil Parikh
  - Dr. Miguel Montero Baker
  - Dr. Bob Beasley